

Product Name: ATG13 Rabbit Monoclonal antibody
Catalog #: AMRe01441



Summary

Production Name	ATG13 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	ATG13
Alternative Names	Autophagy-related protein 13
Gene ID	9776
SwissProt ID	O75143.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 57 kDa; Observed MW: 72 kDa

Background

Product Name: ATG13 Rabbit Monoclonal antibody
Catalog #: AMRe01441




Autophagy factor required for autophagosome formation and mitophagy.


Research Area

Cell Biology

Image Data



Western blot analysis of ATG13 in Raw264.7 lysates using ATG13 antibody.



Western blot analysis of ATG13 in K562, C6, HeLa lysates using ATG13 antibody.

Note

For research use only.