

Product Name: CENPA Rabbit Monoclonal Antibody
Catalog #: AMRe01442



Summary

Production Name	CENPA Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	CENPA
Alternative Names	Centromere autoantigen A; Centromere protein A; CENP-A
Gene ID	1058
SwissProt ID	P49450.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Background

Histone H3-like nucleosomal protein that is specifically found in centromeric nucleosomes.


Product Name: CENPA Rabbit Monoclonal Antibody
Catalog #: AMRe01442



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of CENPA in K562, HeLa lysates using CENPA antibody.

Note

For research use only.