

Product Name: Neurofascin Rabbit Monoclonal Antibody**Catalog #: AMRe01464**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.65mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 138 kDa; Observed MW: 186,155,75 kDa

Antigen Information

Gene Name	Nfasc
Alternative Names	Neurofascin
Gene ID	269116.0
SwissProt ID	Q810U3
Immunogen	Recombinant protein of mouse Neurofascin


Background

Cell adhesion, ankyrin-binding protein which may be involved in neurite extension, axonal guidance, synaptogenesis, myelination and neuron-glial cell interactions.


Research Area

Neuroscience

Image Data



Western blot analysis of Neurofascin in Human placenta lysates using Neurofascin antibody.



Western blot analysis of Neurofascin in mouse Brain, rat Brain, COS7, HeLa lysates using Neurofascin antibody.