

Product Name: SNX5 Rabbit Monoclonal Antibody
Catalog #: AMRe01465



Summary

Production Name	SNX5 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	SNX5
Alternative Names	Sorting nexin-5
Gene ID	27131
SwissProt ID	Q9Y5X3.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Background

Product Name: SNX5 Rabbit Monoclonal Antibody
Catalog #: AMRe01465




Involved in several stages of intracellular trafficking.


Research Area

Signal Transduction


Image Data



Immunocytochemistry analysis of SNX5 (green) in Hela using SNX5 antibody, and DAPI (blue).



Western blot analysis of SNX5 in mouse spleen, mouse brain lysates using SNX5 antibody.



Western blot analysis of SNX5 in K562, rat Brain, C6, Hela lysates using SNX5 antibody

Product Name: SNX5 Rabbit Monoclonal Antibody
Catalog #: AMRe01465



Note

For research use only.