

Product Name: ZHX2 Rabbit Monoclonal Antibody
Catalog #: AMRe01478



Summary

Production Name	ZHX2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	ZHX2
Alternative Names	AFP regulator 1; Zinc finger and homeodomain protein 2
Gene ID	22882
SwissProt ID	Q9Y6X8.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 92 kDa; Observed MW: 110 kDa

Background

Product Name: ZHX2 Rabbit Monoclonal Antibody
Catalog #: AMRe01478




Acts as a transcriptional repressor.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ZHX2 in rat Brain, C6, HeLa lysates using ZHX2 antibody.



Immunocytochemistry analysis of ZHX2 (green) in MCF-7 using ZHX2 antibody, and DAPI (blue)

Note

For research use only.