

Product Name: MAGEC2 Rabbit Monoclonal antibody
Catalog #: AMRe01490



Summary

Production Name	MAGEC2 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	MAGEC2
Alternative Names	CT10,MAGE-C2 antigen,MAGE-E1 antigen,MAGEC2,HCA587, MAGEE1
Gene ID	51438
SwissProt ID	Q9UBF1.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 41 kDa; Observed MW: 55 kDa

Background

Product Name: MAGEC2 Rabbit Monoclonal antibody
Catalog #: AMRe01490




Proposed to enhance ubiquitin ligase activity of RING-type zinc finger-containing E3 ubiquitin-protein ligases.


Research Area

-

Image Data



Immunohistochemistry analysis of paraffin-embedded mouse testis using MAGEC2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of MAGEE1 in K562 lysates using MAGEE1 antibody.

Note

For research use only.