

Product Name: TGF beta 1 Rabbit Monoclonal Antibody
Catalog #: AMRe01510



Summary

Production Name	TGF beta 1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	Tgfb1
Alternative Names	TGF beta 1; TGFB; CED; LAP
Gene ID	21803
SwissProt ID	P04202.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 44 kDa

Background

Product Name: TGF beta 1 Rabbit Monoclonal Antibody
Catalog #: AMRe01510




Transforming growth factor beta-1 proprotein: Precursor of the Latency-associated peptide (LAP) and Transforming growth factor beta-1 (TGF-beta-1) chains, which constitute the regulatory and active subunit of TGF-beta-1, respectively.


Research Area

Cardiovascular

Image Data



Western blot analysis of TGF beta 1 in Raw264.7 lysates using TGF beta 1 antibody.



Western blot analysis of TGF beta 1 in Jurkat, C6 lysates using TGF beta 1 antibody.

Note

For research use only.