

Product Name: Glypican 1 Rabbit Monoclonal Antibody**Catalog #: AMRe01512**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.64mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 55 kDa

Antigen Information

Gene Name	GPC1
Alternative Names	glypican
Gene ID	2817
SwissProt ID	P35052
Immunogen	Recombinant protein of human Glypican 1


Background

Cell surface proteoglycan that bears heparan sulfate.


Research Area

Neuroscience

Image Data



Western blot analysis of Glypican 1/GPC1 in Jurkat, Hela lysates using Glypican 1 antibody.



Immunocytochemistry analysis of Glypican 1/GPC1 (green) in Hela using Glypican 1/GPC1 antibody, and DAPI (blue)