

---

**Product Name: Annexin V Rabbit Monoclonal Antibody****Catalog #: AMRe01522**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ICC/IF
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	0.63mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,ICC/IF 1:50-1:200
<b>Molecular Weight</b>	Calculated MW: 36 kDa; Observed MW: 36 kDa

**Antigen Information**

<b>Gene Name</b>	ANXA5
<b>Alternative Names</b>	PP4; ANX5; ENX2; RPRGL3.
<b>Gene ID</b>	308
<b>SwissProt ID</b>	P08758
<b>Immunogen</b>	A synthetic peptide of human Annexin V/ANXA5


**Background**

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.


## Research Area

Cell Biology

## Image Data



Western blot analysis of Annexin V in K562, C6, 3T3, HeLa lysates using Annexin V antibody.



Immunocytochemistry analysis of Annexin V (green) in HeLa using Annexin V antibody, and DAPI (blue)