

Product Name: Rb Rabbit Monoclonal Antibody
Catalog #: AMRe01530



Summary

Production Name	Rb Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	Rb1
Alternative Names	p110-RB1; pRb; pp105
Gene ID	19645.0
SwissProt ID	P13405.

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 105 kDa; Observed MW: 105 kDa

Background

Product Name: Rb Rabbit Monoclonal Antibody
Catalog #: AMRe01530




Tumor suppressor that is a key regulator of the G1/S transition of the cell cycle.

Research Area

Cell Biology

Image Data



Western blot analysis of Rb in C6, Jurkat, HeLa lysates using Rb antibody.

Note

For research use only.