

Product Name: 14-3-3 theta Rabbit Monoclonal Antibody
Catalog #: AMRe01593

Summary

Production Name	14-3-3 theta Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	YWHAQ
Alternative Names	YWHAQ; 14-3-3 protein theta; 14-3-3 protein T-cell; 14-3-3 protein tau; Protein HS1
Gene ID	10971
SwissProt ID	P27348.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 28 kDa; Observed MW: 28 kDa

Product Name: 14-3-3 theta Rabbit Monoclonal Antibody
Catalog #: AMRe01593


Background

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.


Research Area

Neuroscience

Image Data



Western blot analysis of 14-3-3 Theta in rat brain, rat lung lysates using 14-3-3 theta antibody.



Western blot analysis of 14-3-3 Theta/Tau in 3T3, HeLa lysates using 14-3-3 Theta/Tau antibody.

Note

For research use only.