

Summary

Production Name	AAMP Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	AAMP
Alternative Names	Angio-associated migratory cell protein
Gene ID	14
SwissProt ID	Q13685.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Background

Product Name: AAMP Rabbit Monoclonal Antibody
Catalog #: AMRe01597




Plays a role in angiogenesis and cell migration. In smooth muscle cell migration, may act through the RhoA pathway.


Research Area

Cardiovascular

Image Data



Western blot analysis of AAMP in mouse heart lysates using AAMP antibody.



Western blot analysis of AAMP in C6, Hela lysates using AAMP antibody.

Note

For research use only.