

Product Name: ACAT1 Rabbit Monoclonal Antibody
Catalog #: AMRe01601



Summary

Production Name	ACAT1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	ACAT1
Alternative Names	T2; MAT; ACAT; THIL
Gene ID	38
SwissProt ID	P24752.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 45 kDa

Background

Plays a major role in ketone body metabolism.


Product Name: ACAT1 Rabbit Monoclonal Antibody
Catalog #: AMRe01601



Research Area

Tags & Cell Markers

Image Data



Western blot analysis of ACAT1 in K562, HeLa lysates using ACAT1 antibody.

Note

For research use only.