

Product Name: Adrenodoxin Rabbit Monoclonal antibody
Catalog #: AMRe01615

Summary

Production Name	Adrenodoxin Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	FDX1
Alternative Names	ADX; FDX; LOH11CR1D
Gene ID	2230
SwissProt ID	P10109.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 19 kDa; Observed MW: 14 kDa

Product Name: Adrenodoxin Rabbit Monoclonal antibody
Catalog #: AMRe01615


Background

Essential for the synthesis of various steroid hormones (PubMed:20547883, PubMed:21636783). Participates in the reduction of mitochondrial cytochrome P450 for steroidogenesis (PubMed:20547883, PubMed:21636783). Transfers electrons from adrenodoxin reductase to CYP11A1, a cytochrome P450 that catalyzes cholesterol side-chain cleavage (PubMed:20547883, PubMed:21636783). Does not form a ternary complex with adrenodoxin reductase and CYP11A1 but shuttles between the two enzymes to transfer electrons .


Research Area

Signal Transduction

Image Data



Western blot analysis of Adrenodoxin in 3T3 lysates using Adrenodoxin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Adrenodoxin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.