

Product Name: AER61 Rabbit Monoclonal antibody
Catalog #: AMRe01618



Summary

Production Name	AER61 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	EOGT
Alternative Names	Extracellular O-linked N-acetylglucosamine transferase
Gene ID	285203
SwissProt ID	Q5NDL2.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 62 kDa; Observed MW: 62 kDa

Background

EGF domain-specific O-linked N-acetylglucosamine transferase


Product Name: AER61 Rabbit Monoclonal antibody
Catalog #: AMRe01618



Research Area

Cell Biology

Image Data



Western blot analysis of AER61 in C6, HeLa lysates using AER61 antibody.

Note

For research use only.