

Product Name: AHA1 Rabbit Monoclonal Antibody**Catalog #: AMRe01621**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.3mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 42 kDa

Antigen Information

Gene Name	AHSA1
Alternative Names	AHSA1; C14orf3; HSPC322; Activator of 90 kDa heat shock protein ATPase homolog 1; AHA1; p38
Gene ID	10598
SwissProt ID	O95433
Immunogen	A synthetic peptide of human AHA1

Background


Cochaperone that stimulates HSP90 ATPase activity (By similarity). May affect a step in the endoplasmic reticulum to Golgi

trafficking.


Research Area

Signal Transduction

Image Data



Western blot analysis of AHA1 in rat brain, mouse skin lysates using AHA1 antibody.



Western blot analysis of AHA1 in Hela, HepG2 lysates using AHA1 antibody.