

Product Name: Annexin VII Rabbit Monoclonal Antibody
Catalog #: AMRe01652

Summary

Production Name	Annexin VII Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	ANXA7
Alternative Names	SNX; ANX7; SYNEXIN
Gene ID	310
SwissProt ID	P20073.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 53 kDa

Background

Product Name: Annexin VII Rabbit Monoclonal Antibody
Catalog #: AMRe01652




Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.


Research Area

Signal Transduction

Image Data



Western blot analysis of Annexin VII in K562, rat Brain lysates using Annexin VII antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Annexin VII antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.