

---

**Product Name: Annexin VII Rabbit Monoclonal antibody****Catalog #: AMRe01653**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,IP
<b>Reactivity</b>	Human,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal Antibody
<b>Form</b>	Liquid
<b>Concentration</b>	0.18mg/ml. The concentration of this product may be batch-dependent.
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purified

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 53 kDa; Observed MW: 51,47 kDa

**Antigen Information**

<b>Gene Name</b>	ANXA7
<b>Alternative Names</b>	SNX; ANX7; SYNEXIN
<b>Gene ID</b>	310
<b>SwissProt ID</b>	P20073
<b>Immunogen</b>	Recombinant protein of human Annexin VII


**Background**

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.


## Research Area

Signal Transduction

## Image Data



Western blot analysis of Annexin VII in K562, rat Brain, HeLa lysates using Annexin VII antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Annexin VII antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.