

Product Name: ATG12 Rabbit Monoclonal antibody
Catalog #: AMRe01661



Summary

Production Name	ATG12 Rabbit Monoclonal antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	ATG12
Alternative Names	APG12-like; APG12L; ATG12; Autophagy 12; Autophagy-related protein 12; FBR93; HAPG12; Ubiquitin-like protein ATG12
Gene ID	9140
SwissProt ID	O94817.

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 15 kDa; Observed MW: 55 kDa

Product Name: ATG12 Rabbit Monoclonal antibody
Catalog #: AMRe01661




Background

Ubiquitin-like protein required for autophagy. Conjugated to ATG3 and ATG5.

Research Area

Cell Biology

Image Data



Western blot analysis of Apg12 in K562, 3T3, HeLa lysates using ATG12 antibody.

Note

For research use only.