
Product Name: ARHGAP25 Rabbit Monoclonal Antibody**Catalog #: AMRe01674**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 73 kDa; Observed MW: 73 kDa

Antigen Information

Gene Name	ARHGAP25
Alternative Names	Rho GTPase activating protein 25; KAIA0053
Gene ID	9938
SwissProt ID	P42331
Immunogen	Recombinant protein of human ARHGAP25


Background

GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

Research Area

Signal Transduction

Image Data



Western blot analysis of ARHGAP25 in rat brain, Raoms, Jurkat lysates using ARHGAP25 antibody.