

Product Name: C3orf38 Rabbit Monoclonal antibody
Catalog #: AMRe01745



Summary

Production Name	C3orf38 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	C3orf38
Alternative Names	Uncharacterized protein C3orf38
Gene ID	285237
SwissProt ID	Q5JPI3.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Background

May be involved in apoptosis regulation.


Product Name: C3orf38 Rabbit Monoclonal antibody
Catalog #: AMRe01745



Research Area

Cell Biology

Image Data



Western blot analysis of C3orf38 in K562, HeLa lysates using C3orf38 antibody.

Note

For research use only.