

**Product Name: CaMKII delta Rabbit Monoclonal Antibody**  
**Catalog #: AMRe01751**

---

## Summary

|                        |                                         |
|------------------------|-----------------------------------------|
| <b>Production Name</b> | CaMKII delta Rabbit Monoclonal Antibody |
| <b>Description</b>     | Rabbit Monoclonal antibody              |
| <b>Host</b>            | Rabbit                                  |
| <b>Application</b>     | WB                                      |
| <b>Reactivity</b>      | Human,Rat                               |

## Performance

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                         |
| <b>Modification</b> | Unmodified                                                                                           |
| <b>Isotype</b>      | IgG                                                                                                  |
| <b>Clonality</b>    | Monoclonal                                                                                           |
| <b>Form</b>         | Liquid                                                                                               |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.             |
| <b>Buffer</b>       | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein |
| <b>Purification</b> | Affinity Purification                                                                                |

## Immunogen

|                          |                                                                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CAMK2D                                                                                                                                           |
| <b>Alternative Names</b> | CAM kinase 2 delta; CAM kinase II delta; CaM kinase II delta subunit; CaM kinase II subunit delta; CAMK 2d; CaMK-II delta subunit; CAMK2D; CAMKI |
| <b>Gene ID</b>           | 817                                                                                                                                              |
| <b>SwissProt ID</b>      | Q13557.                                                                                                                                          |

## Application

|                         |                                            |
|-------------------------|--------------------------------------------|
| <b>Dilution Ratio</b>   | WB: 1:500-1:1000                           |
| <b>Molecular Weight</b> | Calculated MW: 56 kDa; Observed MW: 56 kDa |

**Product Name: CaMKII delta Rabbit Monoclonal Antibody**  
**Catalog #: AMRe01751**

---


## Background

CaM-kinase II (CAMK2) is a prominent kinase in the central nervous system that may function in long-term potentiation and neurotransmitter release.

## Research Area

Neuroscience

## Image Data



Western blot analysis of CaMKII delta in rat Brain, C6, HeLa lysates using CaMKII delta antibody.

## Note

For research use only.