

Product Name: Cathelicidin Rabbit Monoclonal antibody
Catalog #: AMRe01766

Summary

Production Name	Cathelicidin Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	CAMP
Alternative Names	LL37; CAP18; CRAMP; HSD26; CAP-18; FALL39; FALL-39
Gene ID	820
SwissProt ID	P49913.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 19 kDa; Observed MW: 19 kDa

Background


Binds to bacterial lipopolysaccharides (LPS), has antibacterial activity.

Product Name: Cathelicidin Rabbit Monoclonal antibody
Catalog #: AMRe01766

Research Area

Microbiology

Image Data



Western blot analysis of Cathelicidin/CLP in Human serum lysates using Cathelicidin antibody.

Note

For research use only.