

Summary

Production Name	CD127 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	IL7R
Alternative Names	IL7R; Interleukin-7 receptor subunit alpha; IL-7 receptor subunit alpha; IL-7R subunit alpha; IL-7R-alpha; IL-7RA; CDw127; CD antigen CD127
Gene ID	3575
SwissProt ID	P16871.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 52 kDa; Observed MW: 70 kDa

Background

Product Name: CD127 Rabbit Monoclonal Antibody
Catalog #: AMRe01775




Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

Research Area

Immunology

Image Data



Western blot analysis of CD127 in K562 lysates using CD127 antibody.

Note

For research use only.