

Product Name: CD134/OX40L Receptor Rabbit Monoclonal Antibody**Catalog #: AMRe01776**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,IP
Reactivity	Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1.36mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 30 kDa; Observed MW: 45 kDa

Antigen Information

Gene Name	Tnfrsf4
Alternative Names	OX40 antigen; OX40L receptor
Gene ID	22163.0
SwissProt ID	P47741
Immunogen	Recombinant protein of mouse CD134


Background

Receptor for TNFSF4/OX40L/GP34. Is a costimulatory molecule implicated in long-term T-cell immunity (By similarity).

Research Area

Immunology

Image Data



Western blot analysis of CD134 in mouse spleen lysates using CD134/OX40L Receptor antibody.