

Product Name: CD272 Rabbit Monoclonal Antibody
Catalog #: AMRe01784



Summary

Production Name	CD272 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	Btla
Alternative Names	A630002H24; B-and T-lymphocyte-associated protein
Gene ID	208154.0
SwissProt ID	Q7TSA3.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 55-80 kDa

Background

Product Name: CD272 Rabbit Monoclonal Antibody
Catalog #: AMRe01784




Lymphocyte inhibitory receptor which inhibits lymphocytes during immune response.


Research Area

Immunology


Image Data



Immunohistochemistry analysis of paraffin-embedded mouse spleen using CD272 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD272/BTLA in 3T3 lysates using CD272/BTLA antibody.



Western blot analysis of CD272/BTLA in K562, HeLa lysates using CD272/BTLA antibody

Product Name: CD272 Rabbit Monoclonal Antibody
Catalog #: AMRe01784



Note

For research use only.