

Summary

Production Name	CD3G Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	CD3G
Alternative Names	CD3 gamma; CD3g; T3G
Gene ID	917
SwissProt ID	P09693.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

Background

Product Name: CD3G Rabbit Monoclonal Antibody
Catalog #: AMRe01794




The CD3 complex mediates signal transduction.


Research Area

Immunology

Image Data



Western blot analysis of CD3G in Jurkat lysates using CD3G antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD3G antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.