

Product Name: CDC25 Rabbit Monoclonal Antibody
Catalog #: AMRe01803



Summary

Production Name	CDC25 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	RASGRF1
Alternative Names	GNRP; GRF1; CDC25; GRF55; CDC25L; H-GRF55; PP13187; ras-GRF1
Gene ID	5923
SwissProt ID	Q13972.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 145 kDa; Observed MW: 145 kDa

Background

Promotes the exchange of Ras-bound GDP by GTP.


Product Name: CDC25 Rabbit Monoclonal Antibody
Catalog #: AMRe01803



Research Area

Signal Transduction

Image Data



Western blot analysis of CDC25 in rat Brain lysates using CDC25 antibody.

Note

For research use only.