

Product Name: Chromogranin A Rabbit Monoclonal Antibody
Catalog #: AMRe01824

Summary

Production Name	Chromogranin A Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	CHGA
Alternative Names	beta Granin; Chromogranin A; ChrA; SP-I; Vasostatin I; Vasostatin-2; Vasostatin II; Pancreastatin; CHGA
Gene ID	1113
SwissProt ID	P10645.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 92 kDa

Product Name: Chromogranin A Rabbit Monoclonal Antibody
Catalog #: AMRe01824




Background

Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.

Research Area

Neuroscience

Image Data



Western blot analysis of Chromogranin A in rat Brain lysates using Chromogranin A antibody.

Note

For research use only.