
Product Name: CLPTM1 Rabbit Monoclonal Antibody**Catalog #: AMRe01835**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 76 kDa; Observed MW: 90 kDa

Antigen Information

Gene Name	CLPTM1
Alternative Names	Cleft lip and palate transmembrane protein 1
Gene ID	1209
SwissProt ID	O96005
Immunogen	A synthetic peptide of human CLPTM1


Background

May play a role in T-cell development.

Research Area

Cell Biology

Image Data



Western blot analysis of CLPTM1 in 293T, U87-MG, Jurkat lysates using CLPTM1 antibody.