

Product Name: COG6 Rabbit Monoclonal Antibody
Catalog #: AMRe01838



Summary

Production Name	COG6 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	COG6
Alternative Names	COD2; SHNS; CDG2L
Gene ID	57511
SwissProt ID	Q9Y2V7.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 73 kDa; Observed MW: 73 kDa

Background

Required for normal Golgi function.


Product Name: COG6 Rabbit Monoclonal Antibody
Catalog #: AMRe01838



Research Area

Tags & Cell Markers

Image Data



Western blot analysis of COG6 in K562, 3T3 lysates using COG6 antibody.

Note

For research use only.