

Product Name: Ct1p1 Rabbit Monoclonal Antibody**Catalog #: AMRe01858**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 91 kDa; Observed MW: 100-145 kDa

Antigen Information

Gene Name	BCL11A
Alternative Names	EVI9; CTIP1; ZNF856; HBFQTL5; BCL11A-L; BCL11A-S; BCL11a-M; BCL11A-XL
Gene ID	53335
SwissProt ID	Q9H165
Immunogen	Recombinant protein of human Ct1p1

Background


Functions as a myeloid and B-cell proto-oncogene. May play important roles in leukemogenesis and hematopoiesis. An essential factor in lymphopoiesis, is required for B-cell formation in fetal liver. May function as a modulator of the

transcriptional repression activity of ARP1.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of CtIP1 in Ramos, Raji lysates using CtIP1 antibody.