

Product Name: GSDME Rabbit Monoclonal Antibody
Catalog #: AMRe01905



Summary

Production Name	GSDME Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GSDME
Alternative Names	DFNA5; ICERE-1; GSDME
Gene ID	1687
SwissProt ID	O60443.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 55 kDa; Observed MW: 55 kDa

Background

Plays a role in the TP53-regulated cellular response to DNA damage probably by cooperating with TP53

Product Name: GSDME Rabbit Monoclonal Antibody
Catalog #: AMRe01905




(PubMed:16897187, PubMed:18223688).

Research Area

Neuroscience

Image Data



Western blot analysis of DFNA5/GSDME in U251MG, U2OS, HepG2, HeLa lysates using GSDME antibody.

Note

For research use only.