

**Product Name: DHRS1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe01909**



## Summary

<b>Production Name</b>	DHRS1 Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	DHRS1
<b>Alternative Names</b>	dehydrogenase/reductase 1; SDR19C1
<b>Gene ID</b>	115817
<b>SwissProt ID</b>	Q96LJ7.

## Application

<b>Dilution Ratio</b>	WB: 1:500-1:1000
<b>Molecular Weight</b>	Calculated MW: 34 kDa; Observed MW: 34 kDa

## Background

NADPH-dependent oxidoreductase which catalyzes the reduction of steroids (estrone, androstene-3,17-dione and

**Product Name: DHRS1 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe01909**




cortisone) as well as prostaglandin E1, isatin and xenobiotics in vitro.

## Research Area

Cell Biology

## Image Data



Western blot analysis of DHRS1 in C6 lysates using DHRS1 antibody.

## Note

For research use only.