

Product Name: ENO3 Rabbit Monoclonal Antibody**Catalog #: AMRe01948**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.64mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 47 kDa; Observed MW: 47 kDa

Antigen Information

Gene Name	ENO3
Alternative Names	enolase 3; MSE; GSD13
Gene ID	2027
SwissProt ID	P13929
Immunogen	A synthetic peptide of human ENO3


Background

Appears to have a function in striated muscle development and regeneration.


Research Area

Signal Transduction

Image Data



Western blot analysis of ENO3 in rat Brain, C6, 3T3, Hela lysates using ENO3 antibody.



Immunocytochemistry analysis of ENO1/ENO2/ENO3(green) in Hela using ENO1/ENO2/ENO3 antibody, and DAPI(blue)