

Product Name: FAM50A Rabbit Monoclonal Antibody**Catalog #: AMRe01972**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Human,Mouse,Rat,Hamster
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.45mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Antigen Information

Gene Name	FAM50A
Alternative Names	family with sequence similarity 50 member A; 9F; XAP5; HXC26; HXC-26; DXS9928E
Gene ID	9130
SwissProt ID	Q14320
Immunogen	A synthetic peptide of human FAM50A


Background

May be a DNA-binding protein or transcriptional factor.


Research Area

Cell Biology


Image Data




Western blot analysis of FAM50A/XAP in HeLa, CHO-K1, C6, Jurkat lysates using FAM50A antibody.



Western blot analysis of FAM50A/XAP in Smooth Muscle, HCT116, HeLa lysates using FAM50A/XAP antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using FAM50A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of FAM50A (green) in HeLa using FAM50A antibody, and DAPI (blue).