

Product Name: FCGR1A Rabbit Monoclonal Antibody
Catalog #: AMRe01978



Summary

Production Name	FCGR1A Rabbit Monoclonal Antibody
Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FCGR1A
Alternative Names	Fc-gamma RI; FcRI; Fc-gamma RIA; FcgammaRIa; CD64
Gene ID	2209
SwissProt ID	P12314.

Application

Dilution Ratio	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 43-80 kDa

Background

Product Name: FCGR1A Rabbit Monoclonal Antibody
Catalog #: AMRe01978




High affinity receptor for the Fc region of immunoglobulins gamma. Functions in both innate and adaptive immune responses.

Research Area

Immunology

Image Data



Western blot analysis of FCGR1A in THP-1 lysates using FCGR1A antibody.

Note

For research use only.