

Product Name: Fetuin A Rabbit Monoclonal Antibody
Catalog #: AMRe01981



Summary

Production Name	Fetuin A Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	AHSG
Alternative Names	16G2; A2HS; AHS; Ahsg; FETUA; FETUB; FetuinA; FetuinB; Gugu; HSGA; IRL685; PRO2743
Gene ID	197
SwissProt ID	P02765.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 45-50 kDa

Background

Product Name: Fetuin A Rabbit Monoclonal Antibody
Catalog #: AMRe01981




Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

Research Area

Cardiovascular

Image Data



Western blot analysis of Fetuin A in HepG2, 3T3, U937 lysates using Fetuin A antibody.

Note

For research use only.