

Product Name: FKBP51 Rabbit Monoclonal Antibody**Catalog #: AMRe01988**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Hamster
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 51 kDa

Antigen Information

Gene Name	FKBP5
Alternative Names	P54; AIG6; FKBP51; FKBP54; PPlase; Ptg-10
Gene ID	2289
SwissProt ID	Q13451
Immunogen	A synthetic peptide of human FKBP51


Background

Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

Research Area

Immunology

Image Data



Western blot analysis of FKBP51 in CHO-K1 lysates using FKBP51 antibody.