
Product Name: FTO Rabbit Monoclonal antibody**Catalog #: AMRe02006**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.17mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 58 kDa; Observed MW: 58 kDa

Antigen Information

Gene Name	FTO
Alternative Names	FTO; KIAA1752; Alpha-ketoglutarate-dependent dioxygenase FTO; Fat mass and obesity-associated protein
Gene ID	79068
SwissProt ID	Q9C0B1
Immunogen	A synthetic peptide of human FTO

Background


Dioxygenase that repairs alkylated DNA and RNA by oxidative demethylation. Has highest activity towards single-stranded

RNA containing 3-methyluracil, followed by single-stranded DNA containing 3-methylthymine.


Research Area

Neuroscience

Image Data



Western blot analysis of FTO in K562, HeLa lysates using FTO antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using FTO antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.