

Product Name: FYB Rabbit Monoclonal antibody
Catalog #: AMRe02010



Summary

Production Name	FYB Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	FYB1
Alternative Names	ADAP; Fyb; FYN binding protein; p120/p130; SLAP130
Gene ID	2533
SwissProt ID	O15117.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 85 kDa; Observed MW: 130,120 kDa

Background

Acts as an adapter protein of the FYN and LCP2 signaling cascades in T cells. Modulates the expression of interleukin-2 (IL-

Product Name: FYB Rabbit Monoclonal antibody
Catalog #: AMRe02010




2). Involved in platelet activation. Prevents the degradation of SKAP1 and SKAP2.

Research Area

Signal Transduction

Image Data



Western blot analysis of FYB in Jurkat lysates using FYB antibody.

Note

For research use only.