

Product Name: GLB1 Rabbit Monoclonal Antibody
Catalog #: AMRe02034



Summary

Production Name	GLB1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GLB1
Alternative Names	GLB1; ELNR1; Beta-galactosidase; Acid beta-galactosidase; Lactase; Elastin receptor 1
Gene ID	2720
SwissProt ID	P16278.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 76 kDa; Observed MW: 80 kDa

Background

Cleaves beta-linked terminal galactosyl residues from gangliosides, glycoproteins, and glycosaminoglycans.


Product Name: GLB1 Rabbit Monoclonal Antibody
Catalog #: AMRe02034



Research Area

Tags & Cell Markers

Image Data



Western blot analysis of GLB1 in K562, C6, 3T3, HeLa lysates using GLB1 antibody.

Note

For research use only.