

Product Name: GLI3 Rabbit Monoclonal Antibody**Catalog #: AMRe02036**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000, ICC/IF 1:50-1:200
Molecular Weight	Calculated MW: 170 kDa; Observed MW: 200 kDa

Antigen Information

Gene Name	GLI3
Alternative Names	GLI3-190; GLI3FL; GLI3-83; GLI3 C-terminally truncated form
Gene ID	463369.0
SwissProt ID	Q5IS56
Immunogen	


Background

Has a dual function as a transcriptional activator and a repressor of the sonic hedgehog (Shh) pathway, and plays a role in limb development.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Gli3 in 293 lysates using GLI3 antibody.