

**Product Name: Glutathione Reductase Rabbit
Monoclonal Antibody**
Catalog #: AMRe02043

Summary

Production Name	Glutathione Reductase Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GSR
Alternative Names	GR; HEL-75; HEL-S-122m
Gene ID	2936
SwissProt ID	P00390.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 56 kDa; Observed MW: 51 kDa

Background


**Product Name: Glutathione Reductase Rabbit
Monoclonal Antibody
Catalog #: AMRe02043**

Maintains high levels of reduced glutathione in the cytosol.

Research Area

Signal Transduction

Image Data



Western blot analysis of Glutathione Reductase in HeLa, Jurkat lysates using Glutathione Reductase antibody.

Note

For research use only.