

**Product Name: Phospho-Glycogen synthase (Ser641)
Rabbit Monoclonal Antibody
Catalog #: AMRe02049**

Summary

Production Name	Phospho-Glycogen synthase (Ser641) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GYS1
Alternative Names	GYS1; GYS; Glycogen [starch] synthase; muscle
Gene ID	2997
SwissProt ID	P13807.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 84 kDa; Observed MW: 84 kDa

**Product Name: Phospho-Glycogen synthase (Ser641)
Rabbit Monoclonal Antibody
Catalog #: AMRe02049**




Background

Transfers the glycosyl residue from UDP-Glc to the non-reducing end of alpha-1,4-glucan. Allosteric activation by glucose-6-phosphate. Phosphorylation reduces the activity towards UDP-glucose. When in the non-phosphorylated state, glycogen synthase does not require glucose-6-phosphate as an allosteric activator; when phosphorylated it does.


Research Area

Signal Transduction

Image Data



Western blot analysis of Glycogen synthase 1/GYS1 (Phospho-Ser641) in C6, 3T3 lysates using Phospho-Glycogen synthase (Ser641) antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Phospho-Glycogen synthase (Ser641) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.