

Product Name: GlyT2 Rabbit Monoclonal antibody**Catalog #: AMRe02050**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 87 kDa; Observed MW: 87 kDa

Antigen Information

Gene Name	SLC6A5
Alternative Names	NET1; GLYT2; HKPX3; GLYT-2
Gene ID	9152
SwissProt ID	Q9Y345
Immunogen	A synthetic peptide of human Glyt2


Background

Sodium- and chloride-dependent glycine transporter 2

Research Area

Neuroscience

Image Data



Western blot analysis of Glyt2 in Jurkat, C6 lysates using GlyT2 antibody.