

Product Name: GOLGA1 Rabbit Monoclonal Antibody
Catalog #: AMRe02055



Summary

Production Name	GOLGA1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GOLGA1
Alternative Names	Golgin subfamily A member 1 (Golgin-97)
Gene ID	2800
SwissProt ID	Q92805.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 88 kDa

Background

The Golgi-associated protein golgin A1 (GOLGA1, golgin-97) was first isolated as a Golgi complex autoantigen associated

Product Name: GOLGA1 Rabbit Monoclonal Antibody
Catalog #: AMRe02055




with the autoimmune disorder Sjogren's syndrome. The golgin-97 protein contains a carboxy-terminal GRIP domain and is a commonly used trans-Golgi network (TGN) marker.

Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Golgin 97 in HeLa lysates using GOLGA1 antibody.

Note

For research use only.