

Product Name: GST3 Rabbit Monoclonal Antibody
Catalog #: AMRe02066



Summary

Production Name	GST3 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	GSTP1
Alternative Names	GSTP1; FAEES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1
Gene ID	2950
SwissProt ID	P09211.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Product Name: GST3 Rabbit Monoclonal Antibody
Catalog #: AMRe02066




Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.


Research Area

Cardiovascular

Image Data



Western blot analysis of GST3/GST pi in HepG2, Hela lysates using GST3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using GST3/GST pi antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.