
Product Name: HBEGF Rabbit Monoclonal Antibody**Catalog #: AMRe02071**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Antigen Information

Gene Name	HBEGF
Alternative Names	DTR; DTS; DTSF; HB-EGF; HBEGF; HEGFL
Gene ID	1839
SwissProt ID	Q99075
Immunogen	


Background

May be involved in macrophage-mediated cellular proliferation. It is mitogenic for fibroblasts and smooth muscle but not endothelial cells.

Research Area

Signal Transduction

Image Data



Western blot analysis of HBEGF in K562, C6, 3T3, HeLa lysates using HBEGF antibody.