

---

**Product Name: HDAC3 Rabbit Monoclonal Antibody****Catalog #: AMRe02076**

For research use only.

**Summary**

<b>Description</b>	Recombinant rabbit monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,IP
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
<b>Molecular Weight</b>	Calculated MW: 49 kDa; Observed MW: 49 kDa

**Antigen Information**

<b>Gene Name</b>	HDAC3
<b>Alternative Names</b>	HDAC3; Histone deacetylase 3; HD3; RPD3-2; SMAP45
<b>Gene ID</b>	8841
<b>SwissProt ID</b>	O15379
<b>Immunogen</b>	A synthetic peptide of human HDAC3

**Background**


HDAC3 is a nuclear and cytoplasmic protein that deacetylates both histone (H2A, H3, H4) and non-histone substrates (RelA, SRY, p53, MEF2, PCAF and p300/CBP). HDAC3 deacetylase activity is stimulated by interactions with the N-CoR and SMRT co-

repressor proteins. Together, these three proteins form a functional complex that represses transcription associated with nuclear hormone receptors and other transcription factors, including Rev-Erb, COUP-TF, DAX1, MAD and Pit-1.

## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of HDAC3 in 3T3, HeLa, rat Brain lysates using HDAC3 antibody.